



CV NOTES

CYANIDE



REACTION I



SAMPLE: HCN
 MEDIUM: 0.1M NaOH
 CONC.: 0.1 mg/mL
 RATE: 200 mV/sec
 ELECTRODE: Silver Planar
 REFERENCE: RE-1, Ag/AgCl/3M KCl
 MODEL: CV-1A

Accurate detection of cyanide is important in many ways. For example, cyanide is an essential ingredient in electrochemical plating processes. It is also a highly toxic anion, causing inactivation of certain respiratory enzymes.

Cyanide can be detected by a novel electrochemical technique utilizing a silver electrode. It is thought that metallic silver is oxidized from the electrode surface to form a cyano complex via Reaction I. Detection by this method appears to be extremely sensitive and readily reversible. Other anions, such as the halides and thiocyanide, are detectable by this method as well.



2701 Kent Ave
 West Lafayette
 Indiana 47906