

## Histamine in Rat Microdialysates by Liquid Chromatography / Electrochemical Detection

1011

### Purpose

Determination of histamine in microdialysates of rat anterior hypothalamus.




Figure 1. Structure of histamine.

Histamine (F1) is a heterocyclic primary amine known for its mediation of allergic response [1]. Stored in mast cells and peripheral blood basophils, its release plays a role in allergy, inflammation, gastric acid secretion, microcirculation and neurotransmission [2].

### Existing Methods

Radioenzymatic assays, fluorometry, and LC with fluorescence or electrochemical detection (both require derivatization). The method below, modified from [3], relies on derivatization with o-phthalaldehyde (OPA).

### Conditions

**System:** [BASi Electrochemical Detector](#) Package with a HPLC pump configured for microbore chromatography.

**Electrochemical Detector Electrode:** Radial-flow 3 mm glassy carbon ([MF-1095](#))

**Potential:** +700 mV vs. Ag/AgCl

**Column:** 3 x 100 mm C<sub>18</sub> (MF-8954)

**Column Temperature:** 35 °C

**Mobile Phase:** 5.9 g NaH<sub>2</sub>PO<sub>4</sub> H<sub>2</sub>O, 1.1 g Na<sub>2</sub>HPO<sub>4</sub>, 186 mg EDTA brought to 500 mL H<sub>2</sub>O, adjusted to pH 6.4. Combine 480 mL with 240 mL acetonitrile and 280 mL methanol.

**Flow Rate:** 1 mL/min

**Loop Size:** 20 µL

**Injection Volume:** 10 µL

### Notes

Samples were derivatized with OPA using Reagent A of the [BASi Amino Acid Kit](#) (MF-8905), using 20 µL sample and 4 µL Reagent A, with a 3-minute reaction time.

Separation of histamine in a rat microdialysis sample is shown in F2. A calibration curve for histamine standards in water is shown in F3. The limit of quantitation for real-world microdialysis samples was approximately 15 nM (0.15 pmoles injected).

### References

1. D. Egger, G. Reisbach, L. Hultner, *J. Chromatogr. B* 662 (1994) 103-107.
2. C.M.C.J. van Haaster, W. Engels, P.J.M.R. Lemmens, G. Hornstra, G.J. van der Vusse, *J. Chromatogr. B* 617 (1993) 233-240.
3. T.B. Jensen, P.D. Marley, *J. Chromatogr. B* 670 (1995) 199-207.




Figure 2. Separation of histamine (H) from rat anterior hypothalamus microdialysate. Peak corresponds to 50 nM histamine.




Figure 3. Linearity of histamine standards. R<sup>2</sup> = 0.997