


# CV NOTES

## THYRONINE



SAMPLE: L-tyrosine

MEDIUM: 1 M  $\text{H}_2\text{SO}_4$

CONC: 5 mg/25 ml

RATE: 200 mV/s

ETRODE: CPE

REF: RE-1, Ag/AgCl

MODEL: CV-1A

A:  $2e^-$  oxidation of thyronine producing p-benzoquinone and tyrosine

B: oxidation of the tyrosine generated from thyronine

C: reduction of p-benzoquinone to p-hydroquinone

D: Oxidation of p-hydroquinone produced via wave C (not present on initial positive scan)

Note that waves C and D will appear after traversing wave A. It is not necessary to pass through wave B to generate p-benzoquinone. This can be proven by reducing the positive switching potential to 0.95 V. Peaks B, C, and D correspond with those for authentic samples of tyrosine, benzoquinone, and hydroquinone respectively.



2701 Kent Ave  
West Lafayette  
Indiana 47906